



Disinfectant cleanser concentrate liquid for instruments and surfaces

Ideal for instruments, large surfaces, floors in clinical and healthcare areas

ICEPUR is a concentrated disinfectant cleanser developed for demanding clinical and healthcare environments. It is suitable for regular use across instruments, larger surfaces and floors in areas where reliable hygiene and efficient routines are essential.



ADVANTAGES

ICEPUR combines broad microbiological efficacy with compatibility across many sensitive materials commonly found in healthcare settings. Its concentrated format supports flexible use across different applications, while low VOCs and a non-flammable formulation make it well suited to sensitive professional environments.

- Alcohol-, formaldehyde-, and phenol-free
- MEA/SSMA tested
- Compatible with many sensitive surfaces like artificial leather, aluminium 6061, Polycarbonate (PC), Polymethyl methacrylate (PMMA), Polypropylene (PP), Silicone, Polyvinyl chloride (PVC), EVA foam, Stainless steel 316
- Proven microbiological efficacy with a broad spectrum of activity
- Effective under dirty conditions
- Low VOCs
- Non-flammable
- Fragrance-free (optional with lemon scent)
- Cleans and disinfects in one step



EJ Solutions SA
Via Santa Maria 68a
CH-6596 Gordola

+41 91 743 6761
info@ejsolutions.ch
ejsolutions.ch

**Protect
what matters.**

AREAS OF APPLICATION

Some examples are:

- **Medical care providers:** Institutions offering direct medical care, diagnosis, and treatment.
- **Dental care providers:** Clinics and professionals providing oral and dental healthcare, including preventive and restorative treatments.
- **Specialized & preventive care:** Centers focused on specific fields, rehabilitation, or preventive health.
- **Pharma & medical supply chain:** Producers and distributors of medicines, equipment, and chemicals.
- **Health support & wellness:** Services supporting emergency response, education, and well-being.

MICROBIOLOGICAL EFFICACY

Efficacy	Microbiology	Contact time
Bactericidal EN 13727, EN 16615, EN 14561	Staphylococcus aureus Enterococcus hirae Escherichia coli Proteus mirabilis Pseudomonas aeruginosa	5 min
Yeasticidal EN 13624, EN 16615, EN 14562	Candida albicans	15 min
Virucidal against EN 14476, EN 17111	Hepatitis C	5 min
Limited spectrum virucidal EN 14476, EN 17111	Norovirus Adenovirus	15 min

PRODUCT ITEM LIST

ICEPUR (Concentrate)		
Size	Form	Art. No.
2 l	Bottle	IC_LB02000_HC_EH
5 l	Canister	IC_LC05000_HC_EH
10 l	Canister	IC_LC10000_HC_EH
2 l	Bottle	IC_LB02000_HC_EH/LM
5 l	Canister	IC_LC05000_HC_EH/LM
10 l	Canister	IC_LC10000_HC_EH/LM

LM = with lemon fragrance

Please note that the availability of products may vary in different countries and regions. Contact your local distributor partner for more information. The recommendations regarding our preparation are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility or microbiological efficacy are possible, only in separate and individual cases. Our recommendations are not binding and do not constitute a guarantee. They do not preclude a company's own testing for the intended purpose and process. In this respect we cannot accept any liability. This is in accordance with our general condition of sale and supply.

DIRECTION FOR USE

Dilute to a 5% solution (50 ml ICEPUR per 1 l water). Apply with a cloth/floor mop to non-absorbent, non-porous surfaces or instruments. Keep the entire surface wet for the required contact time, then rinse with water. Surfaces and instruments intended to come into contact with food, feed or fodder must be rinsed several times with potable water after treatment. Instruments may also be treated by immersion. Disconnect electrical equipment before application. Do not treat food, feed or fodder directly. For professional use. Always read the label before use. Report any serious incidents related to the product to the manufacturer and the competent authority.

STORAGE

ICEPUR disinfectant cleanser must be stored in its original, intact packaging, carefully closed after use, in a dry, cool place at a temperature between 15 °C and 30 °C, and protected from direct sunlight.

REGISTRATION

- BAG CHZN7840
- BAuA N-120975
- UFI: DFDF-SV50-ETK5-DMJ1

PRODUCT TYPES

- PT2 (Surface Disinfectants)
- PT4 (Food and Feed Area)

ACTIVE INGREDIENTS

Alkyl (C12-18) dimethylbenzyl ammonium chloride (ADBAC (C12-18)) - 13.8%.

WARNING AND PRECAUTIONARY STATEMENTS

Hazard Statements

H314 Causes severe skin burns and eye damage.

H410 Very toxic to aquatic life with long lasting effects.

Precautionary Statements

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER/doctor.

P501 DISPOSAL: Dispose of contents / partially emptied and empty containers as hazardous waste.



EJ Solutions SA
Via Santa Maria 68a
CH-6596 Gordola

+41 91 743 6761
info@ejsolutions.ch
ejsolutions.ch

**Protect
what matters.**